

# White Glue

## Light-Cure

Orthodontic Adhesive

## Elegant Aesthetics


White Glue keeps clear white after light-curing, which is similar to tooth color and ceramic brackets. This makes it the best choice for ceramic & sapphire brackets and meets the aesthetic demands of people.



## Strong Bonding Strength

Bonding Strength of adhesive can reach 20MPa after light-curing. Even under poor moisture conditions, it can also reach 18MPa (min. Clinical Effective Bonding Strength is 7MPa), which assures the brackets bonding stability and success rate.

\* The brackets used for bonding test is 3B Ceramic Brackets with groove base.



## High-efficient Treatment

► **Long operating time:** The adhesive will not be significantly cured within 2 minutes under normal ambient light. It ensures that the dentists have sufficient operation time to position brackets.

► **Short light-curing time:** The light-curing time can be shortened to 5 seconds, which saves the chair time greatly to have better treatment result.



## Excellent Experience

► **Etchant:** Gel-type without overflowing, moderate viscosity

► **Primer:** Moisture-tolerant, good fluidity, strong permeability, without irritating odor

► **Adhesive:** Moderate viscosity, fluent extrusion, clear and smooth, easy to remove excess adhesive

## White Glue | Light-Cure Orthodontic Adhesive Kit

- ✂ ×2 Adhesive 4.0g/piece
- ✂ ×1 Etchant 3.0ml/piece
- 🍷 ×1 Primer 3.0ml/piece
- 🪥 ×20 Brush
- 🪡 ×2 Syringe



**Hangzhou Xingchen 3B Dental Instrument & Material Co.,Ltd.**



Add: Building 4, No. 16 Xingyuan Road, DaYuan Town, FuYang District, 311421, Hangzhou, China

Tel: +86-571-87655303

E-mail: hz3b.n1@hz3b.com

Website: www.ortho3b.com